


## pGAL-MF Yeast Expression Vector

### P03304



```

          EcoR I          Hind III Cla I          Sal I   Xho I
          -----
GCT GAA GCT GAA TTC GAT ATC AAG CTT ATC GAT ACC GTC GAC CTC GAG
A  E  A  E  F  D  I  K  L  I  D  T  V  D  L  E
----->
MF alpha reading frame

```

#### Vector features

Position	Feature	Position	Feature
Start: 22 End: 453	GAL-L promoter	Start: 3115 End: 4462	2micron origin of replication
Start: 471 End: 737	Mating factor alpha leader sequence	Start: 4580 End: 5440	AmpR resistance gene
Start: 772 End: 1033	CYC1 terminator	Start: 5588 End: 6255	pBluescript origin of replication
Start: 1818 End: 2621	URA3 auxotrophic marker		